


https://github.com/homurachan/SIRM_RELION

1、Triggering SIRM correction during 3D auto-refine

File Jobs Schedules

I/O Reference CTF Optimisation Auto-sampling Helix Compute Running

Import
Motion correction
CTF estimation
Manual picking
Auto-picking
Particle extraction
Subset selection
2D classification
3D initial model
3D classification
3D auto-refine
3D multi-body
CTF refinement
Bayesian polishing
Mask creation
Join star files
Particle subtraction
Post-processing
Local resolution
Block-based Recons
Reconstruction
External

Mask diameter (A): 180  ?

Mask individual particles with zeros? Yes ?

Use solvent-flattened FSCs? No ?

Do SIRM correction No ?

Threshold of the SIRM SNR weight: 0.6 ?

Histogram threshold: 0.0 ?

Iteration of SIRM: 40 ?

Real space mask: ?

Do MWTCF: No ?

MWTCF value: 0.5 ?

SIRM options

MWTCF options

Schedule Check command Run!

1、Triggering SIRM correction during 3D auto-refine

File Jobs Schedules

I/O Reference CTF Optimisation Auto-sampling Helix Compute Running

Import
Motion correction
CTF estimation
Manual picking
Auto-picking
Particle extraction
Subset selection
2D classification
3D initial model
3D classification
3D auto-refine
3D multi-body
CTF refinement
Bayesian polishing
Mask creation
Join star files
Particle subtraction
Post-processing
Local resolution
Block-based Recons
Reconstruction
External

Mask diameter (A): 180

Mask individual particles with zeros? Yes

Use solvent-flattened FSCs? No

Do SIRM correction Yes

Threshold of the SIRM SNR weight: 0.6

Histogram threshold: 0.0

Iteration of SIRM: 40

Real space mask:

Do MWTCF: Yes

MWTCF value: 0.5

Schedule Check command Run!

Real space
threshold. Using
histogram to
determine it after
each iteration

Fourier SNR weight threshold. It
should be set by trial-and-error.

It's recommended to be less than
60 during refinement
Additional mask

Interpolation number less than this (0.5 here)*
average would not account for resolution.


2、 Validation, Compare the results between RELION and CryoSPARC

1、 Convert the cryoSPARC .cs file into the RELION .star file. Recommend to use the csparc2star.py from pyem package.

I suggest to use `relion_star_handler --i input.star --o output.star` that convert the output star into RELION 3.1ver

2、 Refine utilizing SIRM-RELION

3、 Normally, the reconstructions between cryosparc and relion would have a small rotation difference. Should be corrected first.



Use “Validation_preprocess_read_relion_star_rotate_accordingto_chimera_fitting.py” to correct.

More details: <https://github.com/homurachan/Block-based-reconstruction/wiki/How-to-fit-block-into-original-map>

```
E:\>type J003.txt
0.99970015  0.00208527  0.02439808  -0.09048470
-0.00218999  0.99998850  0.00426617  0.37514684
-0.02438890  -0.00431832  0.99969322  0.64161460

E:\>python ./read_relion_star_rotate_accordingto_chimera_fitting.py
<input star> <chimera txt> <pixel size> <output star> <output INT translation txt>

E:\>python ./read_relion_star_rotate_accordingto_chimera_fitting.py J003.star J003.txt 1.04 J003_fit.star INT.txt
Is relion3.0? = 0
[[ 0.99970015  0.00208527  0.02439808]
 [-0.00218999  0.9999885  0.00426617]
 [-0.0243889  -0.00431832  0.99969322]]
```

4、Validation

We only need to correct the starfile that you want to use for SIRM correction. Let's call it `relion_rotted.star`. The results from cryoSPARC named `cs_refine.star`. Run the 'Validation_compare_cryoSPARC_and_relion.py' script.

```
python Validation_compare_cryoSPARC_and_relion.py relion_rotted.star cs_refine.star "Angular threshold" "n for C-n symmetry" "validation_result.star"
```

The "angular threshold" is the maximum allowed orientation difference in deg. If larger than this threshold, the particle would be recognized as a bad one and move to the other file.

Normally, the recommended value of "angular threshold" is 2~3 for small proteins. It can also be inferred from the reported angle accuracy from RELION. Multiply that value by ~3 is generally good, but you should always try before actually use.

The "n for C-n symmetry" is designed to expand the asymmetric unit into the whole sphere. I found that even applied symmetry, sometimes the starfile from cryoSPARC is the whole sphere, which caused problems for calculation. Please check the orientational distribution graph from cryoSPARC.

4、Validation

```
71
72     for i in range(mline, len(instar_line)):
73         if(instar_line[i].split()):
74             imagename=str(instar_line[i].split()[IMG_index])
75             image_serial=int(imagename.split('@')[0])
76             image_filename_tmp=cryosparc_filename(imagename)
77             image_filename=""
78             cryosparc2relion_serial.append([])
79             cryosparc2relion_imagename.append([])
80             cryosparc2relion_line.append([])
81             for mm in range(0, len(image_filename_tmp)):
82                 image_filename+=image_filename_tmp[mm]
83                 image_filename+="_"
84             cryosparc2relion_serial[len(cryosparc2relion_serial)-1]=image_serial
85             cryosparc2relion_imagename[len(cryosparc2relion_imagename)-1]=image_filename
86             cryosparc2relion_line[len(cryosparc2relion_line)-1]=i
87
88     print("Finish reading part1")
89     orirelion_serial=[]
90     orirelion_imagename=[]
91     orirelion_line=[]
92     for i in range(mline2, len(oristar_line)):
93         if(oristar_line[i].split()):
94             imagename=str(oristar_line[i].split()[ori_IMG_index])
95             image_filename_tmp=cryosparc_filename(imagename)
96             image_serial=int(imagename.split('@')[0])
97
```

It should be noted that these two “cryoSPARC_filename” at Validation_compare_cryoSPARC_and_relion.py need to be assessed before use.

If your relionImageName is cryoSPARC like, which includes random serial at the front of the filename, use the cryoSPARC_filename function. Otherwise use the relion_filename function instead.

4、Validation

The outputs are three files. “validation_result.star” stores the angular changes of each particles. It IS NOT the actual star.

validation_result.star_lessthan_”threshold”_change_line.star is the validation picked results. The orientation information is from the FIRST input star of the script. So if you want to keep the results from cryoSPARC, put that cs_refine.star at the front. However, in our tests, the results from SIRM-RELION are generally better.

The other file: validation_result.star_largerthan_”threshold”_change_line.star stores the dropped particles.

4、 Validation

Summary:

- 1、 Refine the same dataset with cryoSPARC and RELION
- 2、 Fix the rotation difference between two software
- 3、 Set an angular threshold. It removes the particles whose angular changes are larger.

5、Mask-Picking

The Mask-Picking process is normally performed AFTER the Validation process.
Assuming we choose “validation_result.star_lessthan_”threshold”_change_line.star”.

1、 utilize `relion_star_handler` from SIRM-RELION to generate BILD from this star.

`relion_star_handler --i validation_result.star_lessthan_”threshold”_change_line.star --o test.bild --write_bild`. It can be followed by --bild_angle_step and --bild_radius which are shown on the “Write bild options”.


2、 Use the MaskPicking_step1_read_relion_bild_produce_mrc.py to generate 3D mrc from the BILD file.
Require to install numpy and mrcfile.

`python MaskPicking_step1_read_relion_bild_produce_mrc.py test.bild 0.9375 “mrc_radius” test.bild.mrc`


The “mrc_radius” in in pixel and should not be larger than 80, or the next few steps would be painfully slow.

5、Mask-Picking

The BILD



The generated mrc




5、Mask-Picking

It's very clear that our 3D mrc needs smoothness. A very easy way is to lowpass. I recommend to use proc3d from EMAN package.

```
`proc3d test.bild.mrc test.bild_lp.mrc apix=0.5 lp=5`
```

The apix is not to be used later. It only marks that the map is low-pass to 1/5 of the Nyquist frequency.



After low-pass, pick a good threshold and erase the area that you don't want. Generate mask by using `relion_mask_create`.

Next step is to generate an Euler angle table to select particles. In all, all particles whose orientation stay outside of the mask will be removed.

5、Mask-Picking

Generate an Euler angle table.

Please use SIRM-RELION v1.1. The command is:

```
`relion_image_handler --i test_lp_erased_mask.mrc --o nouse.mrc --gen_MP_table --MP_table  
out_anglelist.star --MP_stepsize "angle_step"`
```

The angle_step is the stepsize between sampling point. Normally the default value is enough.

Check the info on the relion_image_handler

```
===== additional subtract options =====  
--optimise_scale_subtract (false) : Optimise scale between maps before subtraction?  
--optimise_bfactor_subtract (0.) : Search range for relative B-factor for subtraction (in A^2)  
--mask_optimise_subtract () : Use only voxels in this mask to optimise scale for subtraction  
===== SIRM parameter. Not recommended for initial refine, except for local refinement. =====  
--SIRM (false) : Perform SIRM reconstruction  
--thres_SIRM (0.3) : SIRM threshold in Fourier space  
--Hist_cutoff (0.0) : cut-off of real space restrain  
--thres_hist (0.0) : SIRM, histogram threshold in real space  
--iter_SIRM (30) : Number of iteration of SIRM  
--SIRM_snr_weight () : SNR weight file in RFLOAT  
--max_radius_SIRM (100.0) : SIRM, histogram threshold in real space  
--gen_MP_table (false) : Read MRC file and generate Mask-Picking table. Use together with --MP_table and --MP_stepsize  
--MP_table () : Output name of Mask-Picking table.  
--MP_stepsize (0.9375) : Stepsize of Mask-picking table
```

5、Mask-Picking

Dropping particles using the Euler angle table.

Please use SIRM-RELION v1.1. The command is:

```
relion_star_handler --i validation_result.star_lessthan_"threshold"_change_line.star --o maskpick_result.star  
--MP_write_particles --MP_table out_anglelist.star --MP_stepsize angle_step"
```

Check the info on the relion_star_handler

```
===== !!!!Special option of Block-based reconstruction =====  
--do_invert_BBR_handeriness (false) : Suppose your handerness was correct.  
===== Write bild options =====  
--write_bild (false) : Generate bild from input star file.  
--bild_angle_step (0.9375) : Angle step of bild  
--bild radius (128.0) : Radius of inner sphere of bild  
===== Mask-Picking, read tables write particles. =====  
--MP_write_particles (false) : Enable Mask-Picking. Use together with --MP_table and --MP_stepsize.  
--MP_table () : Mask-Picking table name  
--MP_stepsize (0.9375) : Stepsize of Mask-picking table
```

5、Mask-Picking

Summery:

- 1、 After validation, convert starfile to BILD and 3D mrc.
- 2、 Low-pass the mrc and erase the area that you don't want to use.
- 3、 Generate angle table and select particles.

Note: Use `relion_star_handler` to handle starfile. Use `relion_image_handler` to handler mrcfile.

6、Final reconstruction

It should be noted that the SNR weight was automatically generated during auto-refine. However, there is no such procedure during reconstruction. So we have to generate a FSC weighting file.

Perform relion_postprocess, then use the script convert_star_2_SIRM_weight.py.

```
`python convert_star_2_SIRM_weight.py postprocess.star weight.star`
```

to generate the FSC weight.

Then use SIRM-RELION GUI.

6、Final reconstruction

The screenshot shows a software interface with a menu on the left and a main panel on the right. The menu includes options like 'Import', 'Motion correction', 'CTF estimation', 'Manual picking', 'Auto-picking', 'Particle extraction', 'Subset selection', '2D classification', '3D initial model', '3D classification', '3D auto-refine', '3D multi-body', 'CTF refinement', 'Bayesian polishing', 'Mask creation', 'Join star files', 'Particle subtraction', 'Post-processing', 'Local resolution', 'Block-based Recons', 'Reconstruction' (highlighted), and 'External'. The main panel has tabs for 'Input', 'SIRM Correction', and 'Running'. The 'Input' tab is active, showing fields for 'Input particles' (with a 'Browse' button), 'Symmetry of particles' (set to 'C1'), 'Maxres in Angstrom' (set to '-1'), 'Padding factor' (set to '2'), 'Random subset' (set to '-1'), and 'Do CTF correction' (set to 'Yes'). A red text overlay reads: 'Reconstruction dialog. Maxres limits the Fourier's distance, -1 means using full space.' At the bottom right are three buttons: 'Schedule', 'Check command', and 'Run!'.

File Jobs Schedules

Input SIRM Correction Running

Input particles: ? Browse

Symmetry of particles: C1 ?

Maxres in Angstrom: -1 ?

Padding factor: 2 ?

Random subset: -1 ?

Do CTF correction Yes ?

Reconstruction dialog. Maxres limits the Fourier's distance, -1 means using full space.

Schedule Check command Run!

From our experience, Separately perform SIRM to the odd/even half could potentially leading to FSC raising problem beyond resolution. As a result, we alternatively to choose reconstructing all the particles, and applying FSC weight, b-factor sharpening and low-pass filter manually.

6、Final reconstruction

File Jobs Schedules

Input SIRM Correction Running

Import
Motion correction
CTF estimation
Manual picking
Auto-picking
Particle extraction
Subset selection
2D classification
3D initial model
3D classification
3D auto-refine
3D multi-body
CTF refinement
Bayesian polishing
Mask creation
Join star files
Particle subtraction
Post-processing
Local resolution
Block-based Recons
Reconstruction
External

Do SIRM correction: Yes ?

Filename of the FSC Weight STAR: Weight.star ?

Threshold of the SIRM SNR weight: 1.0 ?

Histogram threshold: 0.0 ?

Iteration of SIRM: 40 ?

Real space mask: ?

Almost the same as in 3D auto-refine. Except for the "Filename ..." since we do not know the FSC weight. It must be provided for reconstruction

Schedule Check command Run!

6、Final reconstruction

After the reconstruction, use the `relion_image_handler` from SIRM-RELION to apply the filters.

```
===== Post-process like b-factor and FSC weighting. Should be used together with --angpix =====  
--postprocess (false) : Do postprocess for single map.  
--post_fsc () : FSC to be applied.  
--lp (-1) : Low-passed frequency  
--bf (0.) : B-factor to be applied
```

The command is:

```
`relion_image_handler --i reconstct_SIRM_full.mrc --o reconstct_SIRM_full_lp3p6_bfM60.mrc  
--angpix 1.42 --postprocess --post_fsc weight.star --lp 3.6 --bf -60`
```

The `post_fsc` is the exact FSC from the “`convert_star_2_SIRM_weight.py`”.